



RSUD KAJEN KAB. PEKALONGAN

Jl. Raya Karang Sari Karanganyar, Pekalongan, Jawa Tengah

0285-385230

E-mail : rsudkajen@pekalongankab.go.id

HASIL PEMERIKSAAN LABORATORIUM

No.RM : 652614 No.Permintaan Lab : PK202410030070
Nama Pasien : TURAH Tgl.Permintaan : 03-10-2024
JK/Umur : P / 74 Th 3 BI 2 Hr Jam Permintaan : 10:18:26
Alamat : LUMENENG 1/3 PANINGGARAN, -, Tgl. Keluar Hasil : 03-10-2024
No.Periksa : 2024/10/03/000516 Jam Keluar Hasil : 15:27:37
Dokter Pengirim : dr. ISMI CAHYA DHELIMA Kamar : KNG-B7, KENANGA

Pemeriksaan	Hasil	Satuan	Nilai Rujukan	Keterangan
Darah Lengkap 5 Diff				
DARAH LENGKAP :				Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Hemoglobin	13.6	gr/dL	11.7 - 15.5	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Hematokrit	41.9	%	35 - 47	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Lekosit	25.64	$10^3 / \mu\text{L}$	3.6 - 11	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Trombosit	203.0	$10^3 / \mu\text{L}$	150 - 440	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Eritrosit	4.69	$10^6 / \mu\text{L}$	3.8 - 5.2	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
M C V	89.3	fl	80 - 100	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
M C H	29.0	pg	26 - 34	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
M C H C	32.5	g/dL	32 - 36	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome

RDW	12.6	%	11.5 - 14.5	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
PDW	15.8	fL	9.0 - 13.0	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
MPV	9.4	fL	7.9 - 11.1	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
HITUNG JENIS (DIFF) :				Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
- Eosinofil	0	%	2 - 4	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
- Basofil	0	%	0 - 1	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
- Neutrofil Segmen	85	%	25 - 60	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
- Limfosit	9	%	25 - 40	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
- Monosit	6	%	2 - 8	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
- NLR	9.76		1 - 3	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Ureum				
Ureum	48.0	mg/dL	< 48	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Creatinin				
Creatinin	1.1	mg/dL	0.45 - 0.75	Obs Dyspneu ec ADHF profil wet-warm Geriatric syndrome
Gula Darah Stick				
Gula Darah Stick	111			Obs Dyspneu ec ADHF profil

wet-warm
Geriatric
syndrome

Obs Dyspneu
ec ADHF profil
wet-warm
Geriatric
syndrome
Obs Dyspneu
ec ADHF profil
wet-warm
Geriatric
syndrome
Obs Dyspneu
ec ADHF profil
wet-warm
Geriatric
syndrome

Panel Elektrolit				
Natrium	134	mmol/l	135 - 147	
Kalium	4.6	mmol/l	3.5 - 5.0	
Calcium	1.11	mmol/L	1 - 1.3	

Catatan : Jika ada keragu-raguan pemeriksaan, diharapkan segera menghubungi laboratorium.

Tgl.Cetak : 03/10/2024 15.33.12

Penanggung Jawab

Petugas Laboratorium



dr. M.R. RATRI TITISARI S, Sp.PK

KENYO K.P, A.Md